

COMPARISON OF THE SBM PRODUCTS FOR THE ASSESSMENT OF DRY EYE



idra



o.s.a.
PLUS VERSION
(software update)



o.s.a.
BASE VERSION

AUTO-INTERFEROMETRY TEST Automatic evaluation of the lipid layer	INTERFEROMETRY TEST Manual evaluation of the lipid layer	INTERFEROMETRY TEST Manual evaluation of the lipid layer
TEAR MENISCUS-HEIGHT Estimation of the tear lm quantity up to 5 values	TEAR MENISCUS-HEIGHT Estimation of the tear lm quantity up to 5 values	TEAR MENISCUS-HEIGHT Evaluation of the tear lm quantity
AUTO-NIBUT Evaluation of tear lm break-up time. Non-invasive and fully automatic with tear topography and graphic of tear stability	AUTO-NIBUT Evaluation of tear lm break-up time. Non-invasive and fully automatic with tear topography and graphic of tear stability	AUTO-NIBUT Evaluation of tear lm non-invasive break-up time: non-invasive and fully automatic
MEIBOGRAPHY Auto-detection of MGs thanks to infrared LEDs and percentage of loss area	MEIBOGRAPHY Auto-detection of MGs thanks to infrared LEDs and percentage of loss area	MEIBOGRAPHY Auto-detection of MGs thanks to infrared LEDs and percentage of loss area
3D MEIBOGRAPHY The revolutionary introduction of 3D Meibomian Gland imaging. This new imaging system provides strong evidence to support your choice of therapy for your patient.	3D MEIBOGRAPHY The revolutionary introduction of 3D Meibomian Gland imaging. This new imaging system provides strong evidence to support your choice of therapy for your patient.	
EYE BLINKING DETECTION	EYE BLINKING DETECTION	
BUT TEST - STAINING TEST With the use of a yellow lter and a blue LED		
BLEPHARITIS AND CYLINDRICAL DANDRUFF With automatic magnification	BLEPHARITIS AND CYLINDRICAL DANDRUFF With an appropriate additional lens	BLEPHARITIS AND CYLINDRICAL DANDRUFF With an appropriate additional lens
PUPILLOMETRY AND WHITE TO WHITE MEASUREMENT In scotopic, mesopic and photopic light	PUPILLOMETRY AND WHITE TO WHITE MEASUREMENT In scotopic, mesopic and photopic light	PUPILLOMETRY AND WHITE TO WHITE MEASUREMENT In scotopic, mesopic and photopic light
MD. VIGO TREATMENT SUGGESTION Possibility to print a report with suggested diagnosis and treatment	MD. VIGO TREATMENT SUGGESTION Possibility to print a report with suggested diagnosis and treatment	MD. VIGO TREATMENT SUGGESTION Possibility to print a report with suggested diagnosis and treatment
REPORT Different typologies of reports to be printed	REPORT Different typologies of reports to be printed	REPORT Different typologies of reports to be printed
TREATMENT MANAGING	TREATMENT MANAGING	TREATMENT MANAGING
LIFESTYLE QUESTIONNAIRE	LIFESTYLE QUESTIONNAIRE	LIFESTYLE QUESTIONNAIRE
COMPARISON WITH ALL INTERNATIONAL GRADING SCALES (Efron, CCLRU, Jervis)	COMPARISON WITH ALL INTERNATIONAL GRADING SCALES (Efron, CCLRU, Jervis)	COMPARISON WITH ALL INTERNATIONAL GRADING SCALES (Efron, CCLRU, Jervis)